

Polypropylene Lab Storage Cabinet



TopAir's high quality lab storage cabinet combines an ergonomic design and premium materials.


The cabinet is made of white polypropylene featuring a high level of corrosive resistance.

The cabinet complies with international standards, protecting lab staff from inhaling harmful chemicals and providing convenient solution, and storage for bottles and cans.

Hood configuration options

- Lab Storage cabinet with fuming duct connection - can be connected to an existing fuming system.
- Lab Storage cabinet with built-in fan & fuming duct connection - an independent unit that channels the airflow outside the building using flexible ducting.
- Lab Storage cabinet with fan & filter system - an independent unit that provides fuming for the Hood interior with no need for ducting/pipes

- Polypropylene structure with high chemical resistance
- Observation windows from hermetically tempered glass
- Door locks
- Ventilation openings at the sides and top of the Hood
- 3 shelves inside the cabinet
- Optional suction fan



Models

CUSTOM SIZES AVAILABLE!

Spec/Model	LFC-AFF-900-PP	LFC-AFF-1200-PP	LFC-AFF-1600-PP
Description	Lab storage cabinet with independent fume filtering systems	Lab storage cabinet with independent fume filtering systems	Lab storage cabinet with independent fume filtering systems
External Dimensions (W x D x H)	900 x 450 x 2100 mm 35.43 x 26.97 x 82.7"	1200 x 450 x 2100 mm 47.2 x 26.97 x 82.7"	1600 x 450 x 2100 mm 63 x 26.97 x 82.7"
Hood Material	White polypropylene, 6 mm Tempered Glass		
Power Supply	110/220V 50/60 Hz		
Filter	Charcoal Filter / HEPA Filter		

With Filtration System

(including carbon/HEPA filters and powerful economical consumption fan)

Spec/Model	LFC-PF-900-PP	LFC-PF-1200-PP	LFC-PF-1600-PP
Description	Lab storage cabinet with duct fuming connection	Lab storage cabinet with duct fuming connection	Lab storage cabinet with duct fuming connection
External Dimensions (W x D x H)	900 x 450 x 2100 mm 35.43 x 26.97 x 82.7"	1200 x 450 x 2100 mm 47.2 x 26.97 x 82.7"	1600 x 450 x 2100 mm 63 x 26.97 x 82.7"
Hood Material	White polypropylene, 6 mm Tempered Glass		
Power Supply	110/220V 50/60 Hz		
Filter	Charcoal Filter / HEPA Filter		

Accessories

LFC-SPT
Polypropylene spill tray

Tel : +1-855-6-TOPAIR (+1-855-686-7247) Email: sales@topairsystems.com

Web: www.topairsystems.com | Headquarters USA: 8912 68th Avenue Forest Hills New York 11375, USA