



Additional Sizes for Metal Fume Cabinet

Model Type	External Dimensions (WxDxH)	Package Size	Weight
FH-120-A	1200x800x2350 mm, 47.3x31.50x92.5"	1280x880x1700 mm + 1280x880x950 mm 50.4x34.6x66.9" + 50.4x34.6x37.4"	260 kg 573 lb
FH-150-A	1500x800x2350 mm, 59.0x31.5x92.5"	1580x880x1700 mm, 1580x880x950 mm 62.2x34.6x66.9", 62.2x34.6x37.4"	310 kg 683.4 lb
FH-180-A	1800x800x2350 mm, 70.9x31.5x92.5"	1880x880x1700 mm + 1880x880x950 mm 74x34.6x66.9"+ 74x34.6x37.4"	360 kg 793.6 lb
Custom	900x800x2350 mm, 35.4x31.50x92.5"	980x880x1700 mm + 980x880x950 mm 38.5x34.6x66.9" + 38.5x34.6x37.4"	210 kg 462.9 lb
Custom	1600x800x2350 mm, 63x31.50x92.5"	1680x880x1700 mm + 1680x880x950 mm 66.1x34.6x66.9" + 66.1x34.6x37.4"	330 kg 727.5 lb
Custom	2000x800x2350 mm, 78.7x31.50x92.5"	2080x880x1700 mm + 2080x880x950 mm 81.8x34.6x66.9" + 81.8x34.6x37.4"	410 kg 903.8 lb
Custom	2400x800x2350 mm, 94.5x31.50x92.5"	2480x880x1700 mm + 2480x880x950 mm 97.6x34.6x66.9" + 97.6x34.6x37.4"	460 kg 1014 lb