



# Additional Sizes for Active Polypropylene Fume Cabinet

Model Type	External Dimensions (WxDxH)	Package Size	Weight
<b>FH-120-PP-ACT</b>	1200x850x2350 mm, 47.3x33.4x92.5"	1280x880x1700 mm + 1280x880x950 mm 50.4x34.6x66.9" + 50.4x34.6x37.4"	280 kg 617 lb
<b>FH-150-PP-ACT</b>	1500x850x2350 mm, 59.0x33.4x92.5"	1580x880x1700 mm + 1580x880x950 mm 62.2x34.6x66.9" + 62.2x34.6x37.4"	330 kg 727 lb
<b>FH-180-PP-ACT</b>	1800x850x2350 mm, 70.9x33.4x92.5"	1880x880x1700 mm + 1880x880x950 mm 74x34.6x66.9 + 74x34.6x37.4	380 kg 837.7 lb
<b>Custom</b>	900x850x2350 mm, 35.4x33.4x92.5"	980x880x1700 mm + 980x880x950 mm 38.5x34.6x66.9 + 38.5x34.6x37.4	230 kg 507 lb
<b>Custom</b>	1600x850x2350 mm, 63x33.4x92.5"	1680x880x1700 mm + 1680x880x950 mm 66.1x34.6x66.9" + 66.1x34.6x37.4"	350 kg 771.6 lb
<b>Custom</b>	2000x850x2350 mm, 78.7x33.4x92.5"	2080x880x1700 mm + 2080x880x950 mm 81.8x34.6x66.9" + 81.8x34.6x37.4"	430 kg 947.9 lb
<b>Custom</b>	2400x850x2350 mm, 94.5x33.4x92.5"	2480x880x1700 mm + 2480x880x950 mm 97.6x34.6x66.9" + 97.6x34.6x37.4"	480 kg 1058.2 lb