

# Airflow Alarm



The AirFlow Alarm is an advanced system for ductless fume hoods and fume Hoods, which monitors the airflow performance, and provides visual and audio alerts upon deviations.

The system features an elegant glass panel, enhanced with a mute button that can silence the alarm at the customer's convenience. The Airflow Alarm offers two models: One alerts on a low air velocity level only, and the other alerts on both a high and low airflow level.

The Airflow Alarm ensures that ductless fume hoods operate in a safe and fully functional mode, for the operator and staff.

- Elegant glass panel
- Simple user friendly interface
- Choice of 2 models
- Maintenance, consulting and replacement parts are conveniently available from AAC Control brand

## Specifications

Part	Catalog No.
Accuracy	+ 0.07 m/s
Visual display	Green and red LEDs
Alarm indications	LED and audible alarm
Input power	230/115 VAC, 50/60hz
Mounting	Semi flush

## Models

Part	Catalog No.
Advanced system providing alerts on low/high-level air-flow in the fume hood.	VA-AFA-LH-CI
Basic system providing alerts on low-level airflow in the fume hood.	VA-AFA-L-CI

