

## Carbon Filters

Catalog no.	Description	Size
CF-90-CR	Carbone filter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-CR	Carbone filter for Ductless fume hood size 120 cm	1100x300x70mm
CF-150-CR	Carbone filter for Ductless fume hood size 150 cm	1400x300x70mm
CF-180-CR	Carbone filter for Ductless fume hood size 180 cm	1700x300x70mm

## Acid Filters

Catalog no.	Description	Size
CF-90-AC	Acid filter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-AC	Acid filter for Ductless fume hood size 90 cm	1100x300x70mm
CF-150-AC	Acid filter for Ductless fume hood size 90 cm	1400x300x70mm
CF-180-AC	Acid filter for Ductless fume hood size 90 cm	1700x300x70mm

## Basic Filters

Catalog no.	Description	Size
CF-90-BC	Basis filter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-BC	Basis filter for Ductless fume hood size 120 cm	1100x300x70mm
CF-150-BC	Basis filter for Ductless fume hood size 150 cm	1400x300x70mm
CF-180-BC	Basis filter for Ductless fume hood size 180 cm	1700x300x70mm

## Formaldehydes Filter

Catalog no.	Description	Size
CF-90-FMD	Formaldehydes filter for ductless fume hood 90 cm	800x300x70mm
CF-120-FMD	Formaldehydes filter for ductless fume hood 90 cm	1100x300x70mm
CF-150-FMD	Formaldehydes filter for ductless fume hood 90 cm	1400x300x70mm
CF-180-FMD	Formaldehydes filter for ductless fume hood 90 cm	1700x300x70mm

## Multi Gas Filters (when using both Acids and Basis)

Catalog no.	Description	Size
CF-90-MG	Multi Gas filter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-MG	Multi Gas filter for Ductless fume hood size 120 cm	1100x300x70mm
CF-150-MG	Multi Gas filter for Ductless fume hood size 150 cm	1400x300*70mm
CF-180-MG	Multi Gas filter for Ductless fume hood size 180 cm	1700*300*70mm

## Pre Filters

Catalog no.	Description	Size
CF-90-FR	Prefilter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-FR	Prefilter for Ductless fume hood size 120 cm	1100x300x70mm
CF-150-FR	Prefilter for Ductless fume hood size 150 cm	1400x300x70mm
CF-180-FR	Prefilter for Ductless fume hood size 180 cm	1700x300*70mm

## HEPA Filters

Catalog no.	Description	Size
CF-90-HP	Hepa filter for Ductless fume hood size 90 cm	800x300x70mm
CF-120-HP	Hepa filter for Ductless fume hood size 120 cm	1100x300x70mm
CF-150-HP	Hepa filter for Ductless fume hood size 150 cm	1400x300x70mm
CF-180-HP	Hepa filter for Ductless fume hood size 180 cm	1700x300x70mm

