

Invant-uk

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175

CERTIFICATE & REPORT NO: INV/EN14175/562-B ISSUE DATE: 24th April 2013

Fume Cupboard Type:

chemtop1200 **CF-120-PP ductless bench type External Dimensions:** *Height* = 2060mm (1260mm + 800mm)Width = 1200mm Depth = **700mm** (max) Internal Dimensions: *Height* = **700mm** (max) Width = 1125mm Depth = 620mm (max)**Baffle Gap Dimensions:** 32x75mmx8mm slots Top: Bottom: 116x75mmx8mm slots Fume Cupboard Manufacturer: TopAir USA 3182 Monterey Drive Merrick NY 11566 USA

Fume Cupboard Flow:Volume flow rate:0.310m³/sTest opening width:1125mmTest opening height:500mmFace velocity:0.520m/sPressure drop:N/A

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

Tested and Certified by: Dr A F Bicen Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A.F.B.