

Invant-uk

## **CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175**

CERTIFICATE & REPORT NO: INV/EN14175/560-B ISSUE DATE: 24th April 2013

## **Fume Cupboard Type:**

chemtop900 **CF-090-PP ductless bench type External Dimensions:** *Height* = 2060mm (1260mm + 800mm)Width = 900mm Depth = **700mm** (max) **Internal Dimensions:** *Height* = **700mm** (max) Width =825mm **620mm** (max) *Depth* = **Baffle Gap Dimensions:** 23x75mmx8mm slots Top: Bottom: 83x75mmx8mm slots

Fume Cupboard Manufacturer: TopAir USA 3182 Monterey Drive Merrick NY 11566 USA

Fume Cupboard Flow:Volume flow rate:0.230m³/sTest opening width:825mmTest opening height:500mmFace velocity:0.518m/sPressure drop:N/A

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

**Tested and Certified by:** Dr A F Bicen Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A.F.B.